


GB

Axtair One®Plus

Axtair Automorpho®Plus

Axtair Axensor® AT12/AT15/AT20

Axtair XXL®



CE



WINNCARE FRANCE - Site ASKLÉ SANTÉ  
200 Rue Charles Tellier - Actiparc de Grézan  
30034 Nîmes Cedex 1  
FRANCE


Tél. : +33 (0)4 66 02 15 15  
Fax. : 33 (0)4 66 02 15 00  
Email : [contact@winncare.com](mailto:contact@winncare.com)  
[www.winncare.com](http://www.winncare.com)



MAJ 21/10/2021 GB

# Axtair One®Plus

	VAXT2/ONE-P	VAXT2/ONE/CIC-P
	VAXT1/POMPE/ONEP	
	VAXT2/MA/ONE	VAXT2/MA/ONE/CIC



**CPR**

30 ≤ 150 Kg  
66 ≤ 330 Lb

< 35 dBA

12 cm


8 h

5,7 kg  
12,5 Lb

195 x 87x 17 cm  
76,7x34,2x 6,7 inch

# Axtair Automorpho®Plus

	VAXT4/AUTO-P VAXT4/CIC-P	VAXT4/AUTO100-P VAXT4/CIC100-P	VAXT4/AUTO120-P VAXT4/CIC120-P
	VAXT1/POMPE/AUTO		
	VAXT/MA/AUTOP VAXT/MA/AUTOP/CIC	VAXT/MA/AUTOP100 VAXT/MA/AUTOP100/CIC	VAXT/MA/AUTOP120 VAXT/MA/AUTOP120/CIC



**CPR**

30 ≤ 165 Kg  
66 ≤ 364 Lb

< 35 dBA

12 cm


8 h

8 < 12kg  
17,6 < 26,4Lb

195 x 87(<117) x 17 cm  
76,7x34,2[<46]x 6,7 inch

# Axtair Axensor® AT12

	VAXT6/AUTO	VAXT6/AUTO-S
	VAXT6/POMPE/AUTO	<del>VAXT6/POMPE/AUTO-S</del>
	VAXT6/MA/AUTO	



**CPR**

30 ≤ 165 Kg  
66 ≤ 364 Lb

< 35 dBA


12 cm


8 h



8 kg  
17,6 Lb

195 x 87x 17 cm  
76,7x34,2x 6,7 inch

# Axtair Axensor® AT15

 30 ≤ 180 Kg  
66 ≤ 396 Lb


	VAXT6/MAX	VAXT6/MAX-S
	VAXT6/POMPE/AUTO	<del>VAXT6/POMPE/AUTO-S</del>
	VAXT6/MA/MAX	




 < 35 dBA




 8 h



 10 kg  
22 Lb

 195 x 87 x 20.5 cm  
76.7x34.2x8 inch

# Axtair Axensor® AT20

 30 ≤ 200 Kg  
66 ≤ 440 Lb


	VAXT6/AT20-80	VAXT6/AT20-90	VAXT6/AT20-80-S	VAXT6/AT20-90-S
	VAXT6/POMPE/AUTO		<del>VAXT6/POMPE/AUTO-S</del>	
	VAXT6/MA/AT20-80		VAXT6/MA/AT20-90	




 < 35 dBA




 8 h



 10kg  
22 Lb

 54 x 54.5 x 56 cm  
21,2 x 21,2 x 22 inch

# Axtair XXL®


 100 ≤ 270 Kg  
220 ≤ 595 Lb


	VAXT/XXL120	VAXT/XXL140
	VAXT/POMPE/XXL	
	VAXT/MA/120XXL	VAXT/MA/140XXL



 < 35 dBA

 8 h

 21 < 23 kg  
46 < 50 Lb

 195 x 117 (<137) x 20.5 cm  
76.7x46 (<54) x 8 inch

## 1. INDICATIONS

### Device intended use

This medical device is intended to be used medically in the treatment and prevention of pressure sores.

### Indications

Prevention and support for the treatment of stage 1 to 4 pressure ulcers (according to medical opinion) for patients who may or may not be up during the day and/or have a "moderate to very high" risk of developing a pressure sore, assessed according to a proven scale and based on clinical judgement. (See diagrams on the inside front cover)

### Contraindications

Patient weight Min<Max. Non-stabilized post-traumatic fractures. For use in hyperbaric chambers and on stretchers.

### Target group of patients and users


Hospitalized, institutionalized or home-based adults of over 146 cm in height, with one or more pressure sores and/or at risk of developing pressure sores due to the transitory or definitive alteration of their condition. These people are cared for by health professionals, assisted by carers when necessary.

### Identify adverse side effects



Inform the competent authority if you consider or have reason to believe that the device poses a serious risk or has been tampered with.

## 2. MEDICAL DEVICE COMPOSITION



- 1 Compressor
- 2 Systems for attaching to the sides of the medical bed
- 3 Power supply cable
- 4 Alternating cell mattress, 2 static cells in the head zone 4 static cells in the foot zone
- 5 Individual deflation valve
- 6 Polyurethane foam base
- 7 Pneumatic connector fitted with a plug for rebalancing support pressure when disconnected from the compressor.
- 8 Identification labels, Dirty/clean label,
- 9 User manual,
- 10 Simplified instructions attached to the side of the compressor.

## 3. CLINICAL BENEFIT, PERFORMANCE, MECHANISM OF ACTION

### Device performance characteristics

- > Operating principle: "mechanical" effect based on the alternating inflation of the overlay mattress cells and the pneumatic management of the applied pressure.
- > The inflation level adjustment is automatic and based on the patient's morphology. No external intervention is required.


> "Dynamic" mode: alternating pressures prevent prolonged vascular compression that can lead to tissue hypoxia.

> Low pressure "static" mode: immobilization (orthopaedic, neurological trauma), local pain, withdrawal phases. This mode is not active when the compressor is connected to a cushion.

> "Care" mode: handling, performance of certain medical acts and transfers. This mode is not active when the compressor is connected to a cushion.

### Projected clinical benefits

- > Maintenance of tissue oxygenation in the anatomical areas in contact with the surface of the support by decreasing the pressure applied to the cutaneous and subcutaneous tissue.








## 4. INSTRUCTIONS FOR USE

### User training and qualification

Users trained by individuals approved by the economic operators, notably in safety and non-compliance reporting aspects

### Device Installation

- › Check compressor and (overlay) mattress compatibility (see "technical specifications" table)



- › Assess the risk of patient entrapment in non-moving parts according to IEC 60601-2-52 (the requirements of figures 201.107, 201.108 and Table 201.101 are excluded)

### Cleaning and disinfection

- › Between each patient.
- › Bio-cleaning or steam process.
- › Surface detergent and disinfectant products conforming to requirements of Regulations (UE) n°648/2004 et n°528/2012.
- › Prohibit the high-pressure jet process.
- › Prohibit colouring, industrial degreasers, abrasive and solvent-based products.

### Preventive maintenance

- › Check the device every 2 years of use or after 17500 hours of operation (Indicator: maintenance key led).
- › Contact the manufacturer or distributor in regard to the AIRCARE maintenance solution (training, software, connection kit, revision kit).
- › See Technical Manual (Downloadable from [www.winncare.com](http://www.winncare.com))



## 5. WARNINGS, PRECAUTIONS FOR USE

### Precautions for use

- › Non-stabilized bone and/or muscle injuries in contact with the support.
- › Initial days of post pressure ulcer surgery (skin graft or flap) [Prefer static low pressure mode].
- › Home-based care without the possibility of medical assistance.
- › Bedridden patient weighing more than 135 kg in a semi-seated position > 45°: check the absence of contact with the bed base via a "judged" test by placing the hand palm upwards between the gluteal area and the support. "Comfort" setting can be used to add air.

### Warning

- › If the device alarm LED is flashing or fixed contact your maintenance department as quickly as possible in order to carry out the appropriate troubleshooting.
- › Installation and commissioning according to the EMC information provided by WINNCARE on request.
- › Only use accessories and cables supplied and/or specified by WINNCARE
- › Observe the storage and operating conditions specified by Winncare.
- › Associate the compressor reference with its support: (technical data tables)

### Action required

#### The support alone is not enough to prevent bed sores:

- › Change position at least every 2 to 3 hours
- › Maintain skin hygiene and avoid maceration
- › In case of incontinence, change protection regularly
- › Ensure a sufficient and appropriate diet
- › Drink regularly and in sufficient quantities
- › Avoid unnecessary thickness and foreign objects between the body and support.

### Notify your doctor or nurse

- › Of any abnormal event (fever, pain, redness, or whitening of the support points of your body with the support)
- › If the required measures for use of the medical device cannot be observed.






Indicates the compressor is switched on





Flashes to indicate inflation of the support  
After switch-off, the patient can be installed




CPR (emergency or withdrawal)  
Open position for inflation or deflation  
Closed position during use.



AXENSOR AT20 - XXL



VAXT4 / VAXT6



AXENSOR AT20 - XXL



Indicator light on: keypad automatically locked after 5 minutes of continuous  
4 second press  
Unlocking: continuous 4 second press.



Care mode (static): mode duration limited to 30 minutes  
Idle mode when the compressor is connected to a cushion.



Care mode indicator flashes 5 minutes before the end. Audible signal emitted  
when triggered.  
After 30 minutes: switches to the previously used mode.



Fixed LED: Low priority alarm. Contact the maintenance department



Flashing LED: Medium priority alarm. Remove the patient.  
Contact the maintenance department



Press button: stops sound alarm  
Moderate alarm: reactivation after 3 minutes.



See Technical Manual (Downloadable on [www.winncare.com](http://www.winncare.com))



Caution, read the user manual and (or) the technical manual.



Category II device (Dual insulation)



BF type electrical device (applied to supports)



Complies with the general requirements of Regulation (EU) 2017/745 relating to medical devices



Caution, electrical and electronic equipment subject to selective waste collection



Manufacture date



Serial number



Batch number



Patient weight range



Warning



Wash with water, T° max 90°C, reduced mechanical action, reduced temperature rinsing, reduced spin.



Bleaching possible, chlorination at 5000 ppm allowed.



Ironing excluded.



Dry cleaning not allowed, use of solvent-based stain remover not allowed.



Tumble drying allowed, moderate temperatures (60°C)

## Use



## Storage



# WIN CARE

## FRANCE

4, Le Pas du Château 85670 Saint Paul Mont Penit - FRANCE  
Tél : +33 (0)2 51 98 55 64

## ESPAGNE

Ctra. Masía del Juez 37b 46909 Torrent Valencia - ESPAGNE  
Tél : +34 9 61 56 55 21

## POLOGNE

ul. Lubska 17PL 68-320Jasień - POLOGNE  
Tél : +48 68 371 10 45

## DANEMARK

Tålborgvej 12CDK-4220Korsør - DANEMARK  
Tél : +45 70 27 37 20

## ROYAUME UNI

Unit 15 Gregory Way Stockport, Cheshire, SK5 7ST - UNITED KINGDOM  
Tél : +44 161 477 7900



Distributor

