



# ALOVA™ CONTOURED



## KEY FEATURES

- ▶ Molded visco elastic foam cushion: density  $\geq 75 \text{ Kg/m}^3$
- ▶ Minimum thickness  $> 6 \text{ cm}$  and integrated coccyx cut-out for a reinforced prevention of ischial and sacrococcygeal pressure ulcer risks
- ▶ External rotator to avoid hip abduction
- ▶ Tilting position to reduce the risk of slipping forward
- ▶ Removable and washable Promust PU integral cover (polyurethane enduction on jersey) with non slip-base and gripping handle
- ▶ 16 sizes and 3 heights available



Mapping of pressure distribution values for 80 kg/172 cm patient.



## INDICATIONS

HELP IN PREVENTING PRESSURE ULCERS FOR PATIENT AT MEDIUM TO HIGH RISK OF PRESSURE ULCER DEVELOPMENT

PATIENT WEIGHT: 20 TO 200 KG/ 3.1 TO 31.4 ST DEPENDING ON MODELS

- ▶ Reinforced prevention of pressure ulcer risks in ischial and sacral area
- ▶ External rotator to ensure proper alignment of the lower limbs
- ▶ Tilting position to reduce the risk of slipping forward
- ▶ Removable, stretchable and breathable cover for a thermoregulation when prolonged sitting
- ▶ Inner cover with non-slip base to ensure the cushion keeps in place
- ▶ Various sizes available to meet different needs

DESCRIPTION	ITEM REFERENCE	WIDTH CM/INCH	DEPTH CM/INCH	MIN/MAX HEIGHT CM/INCH*	PATIENT WEIGHT KG/ST	COVER	FIRE STANDARDS
ALOVA® CURVED STANDARD CUSHION	VCAG3636	36/14.1	36/14.1	6<9- 2.3 < 3.5	20 to 80 Kg/3.1 to 12.5 St	Integral cover in Promust PU (polyurethane coated jersey)  black on top and silver on the sides.  Lower side non-slip in PU/PVC fitted with a carrying handle.	Compliant with standard EN 597-1 and EN 597-2.
	VCAG4343	43/16.9	43/16.9	6<9- 2.3 < 3.5	30 to 100 Kg/4.7 to 15.7 St		
	VCAG4340	43/16.9	40/15.7	6<9- 2.3 < 3.5	30 to 90 Kg/4.7 to 14.1 St		
	VCAG4540	45/17.7	40/15.7	6<9- 2.3 < 3.5	30 to 120 Kg/ 4.7 to 18.8 St		
	VCAG3841	38/14.9	41/16.1	6<9- 2.3 < 3.5	30 to 90 Kg/4.7 to 14.1 St		
	VCAG4041	45/17.7	41/16.1	6<9- 2.3 < 3.5	30 to 90 Kg/4.7 to 14.1 St		
	VCAG4341	43/16.9	41/16.1	6<9- 2.3 < 3.5	30 to 100 Kg/4.7 to 15.7 St		
	VCAG4043	450/19.6	43/16.9	6<9- 2.3 < 3.5	30 to 90 Kg/4.7 to 14.1 St		
	VCAG4544	45/17.7	44/17.3	6<9- 2.3 < 3.5	30 to 120 Kg/ 4.7 to 18.8 St		
	VCAG5044	50/19.6	44/17.3	6<9- 2.3 < 3.5	30 to 135 Kg/4.7 to 21.2 St		
	VCAG4347	43/16.9	47/18.5	6<9- 2.3 < 3.5	30 to 120 Kg/ 4.7 to 18.8 St		
	VCAG5250	52/20.4	50/19.6	7<10 - 2.7<3.9	120 to 160 Kg/18.8 to 25.1 St		
	VCAG5450	54/21.2	50/19.6	7<10 - 2.7<3.9	120 to 160 Kg/18.8 to 25.1 St		
VCAG5650	56/22	50/19.6	7<10 - 2.7<3.9	120 to 200 Kg/18.8 to 31.4 St			
VCAG6050	60/23.6	50/19.6	7<10 - 2.7<3.9	120 to 200 Kg/18.8 to 31.4 St			
ALOVA® CURVED TECHNICAL CUSHION WITH PELVIC SUPPORT	VCAG4343XL	43/16.9	43/16.9	9 < 12 - 3.5<4.7	30 to 110 Kg/4.7 to 17.3 St		
	VGAG5044XL	50/19.6	44/17.3	9 < 12 - 3.5<4.7	30 to 135 Kg/4.7 to 21.2 St		
	VCAG3841BP	38/14.9	41/16.1	6<9- 2.3 < 3.5	30 to 90 Kg/4.7 to 14.1 St		
	VCAG4041BP	45/17.7	41/16.1	6<9- 2.3 < 3.5	30 to 90 Kg/4.7 to 14.1 St		
	VCAG4341BP	43/16.9	41/16.1	6<9- 2.3 < 3.5	30 to 100 Kg/4.7 to 15.7 St		
	VCAG4343BP	43/16.9	43/16.9	6<9- 2.3 < 3.5	30 to 100 Kg/4.7 to 15.7 St		
	VCAG4544BP	45/17.7	44/17.3	6<9- 2.3 < 3.5	30 to 120 Kg/ 4.7 to 18.8 St		
	VCAG5044BP	50/19.6	44/17.3	6<9- 2.3 < 3.5	30 to 135 Kg/4.7 to 21.2 St		
	VCAG4347BP	43/16.9	47/18.5	6<9- 2.3 < 3.5	30 to 120 Kg/ 4.7 to 18.8 St		

\* min. height = cushion height measured at the curve base / max. height =total cushion height  
Cover compatible with detergent-disinfectant products for hospital use.

3 years warranty

Cleaning instructions: Promust PU cover



Photos not contractually binding

FP000044/01-2018



WINNOCARE Group - 200 rue Charles Teulier - Actiparc de Grézan - 30034 Nîmes Cedex 1 France  
Ph : +33 (0)4 66 02 15 15 - Fax : +33 (0)4 66 02 15 00 - Email : contact.askle@winnocare.fr - www.winnocare.fr  
In the interest of constantly improving its products, the WINNOCARE Group reserves the rights to modify the technical specifications of any product without prior notice  
Class I medical device in accordance with European Directive 93/42/EEC

