

CARMEN ROOM



UNDERSTATED AND
RELAXING DESIGN...



CARMEN Atmosphere
Abélia II boards
Cream surround, Light Dragon Tree interior finish

FURNITURE FOR THE CARMEN ROOM



BEDSIDE CABINET

- DIMENSIONS:
W445 x D447 x H681 mm
W538 x D560 x H814 mm with 30 L refrigerator
- ▶ 2 rear mobility castors
 - ▶ 1 drawer with removable tray
 - ▶ 1 shelf with bottle space
 - ▶ Hinge door hung to the left as standard, can be hung to the right as an option



DRESSER

- DIMENSIONS:
W875 x D548 x H758 mm
- ▶ 3 tiroirs sur coulisse métal avec butée anti-chute



SINGLE-DOOR WARDROBE

- DIMENSIONS:
W675 x D548 x H1793 mm
- ▶ Left-hand hinge
 - ▶ 1 shelf
 - ▶ 1 clothes rail



TABLE

- DIMENSIONS:
W875 x D548 x H758 mm



TABLE WITH DRESSER

- DIMENSIONS:
W1416 x D548 x H758 mm
- ▶ One leg for open access, bridge chair access 600 mm wide and 720 mm high
 - ▶ 3 right-hand drawers



TWO-DOOR WARDROBE

- DIMENSIONS:
W1125 x D548 x H1793 mm
- ▶ 1 full-width top shelf
 - ▶ 3 half-shelves
 - ▶ 1 right-hand clothes rail



DESK

- DIMENSIONS:
W1075 x D548 x H758 mm
- ▶ 1 right-hand drawer
 - ▶ Rear access for the bridge chair



CHAIR

- DIMENSIONS:
W470 x D530 x H815 mm
- SITTING:
W470 x D460 x H460 mm
- ▶ Choice of PVC coating



BRIDGE

- DIMENSIONS:
L560 x P560 x H815 mm
- ASSISE:
L440 x P460 x H460 mm
- ▶ Choice of PVC coating
 - ▶ Arm rest height: 68 cm

CARMEN	FURNITURE AND BOARD COLORS		SURROUND, EDGING AND HAND RAIL COLORS	MATERIAL
	STANDARD COLORS	OTHER COLORS AS AN OPTION		
ABELIA II BOARDS	Inside finish Light Beech, surround and handrail rail 1001	Choice of melamine finishes from our color chart, 19 mm thick	Painted steel hand rail: Ral ral 1001 and polyethylene ends Ral ral 1001 Edging: Ral 1001 glued ABS	Painted steel hand rail Chipboard boards 19 mm Melamine-faced ABS edging
CARMEN II BOARDS	Light Beech inside finish, matching edging and footboard: smooth epoxy hand rail and polyethylene ends in RAL 8011 brown	Choice of melamine finishes from our color chart, 30 mm thick	Footboard: Painted steel hand rail: Ral ral 1001 and Polyethylene ends Ral ral 1001 Glued ABS edging, Ral ral 1001 color	Painted steel hand rail 30 mm melamine-faced chipboard boards ABS edging
CARMEN II BOARDS spaces for wooden side rails				
FURNITURE	Light Beech inside finish, ral 1001 surround	Choice of melamine finishes from our color chart, 19 mm thick	Ral 1001 standard color edging	Chipboard unit with 19 mm melamine-face / 1 mm ABS edging Chipboard boards with 19 mm melamine-face / 1 mm ABS edging

1-year warranty*. *Under normal conditions of use and excluding wearing parts - please see instructions of use

Photos not contractually binding

WINCARE Group - 4, Le Pas du Château - 85670 Saint-Paul-Mont-Penit - France
Tél : +33 (0)2 51 98 55 64 - Fax : +33 (0)2 51 98 59 07 - E-mail : info@medicatlantique.fr - www.winnicare.fr

In the interest of constantly improving its products, the WINCARE Group reserves the rights to modify the technical specifications of any product without prior notice

Class I medical device in accordance with European Directive 93/42/EEC

