



## LOW BACK SLING

Supportive half sling, allowing good contact between the user and the carer.

### USAGE

- For users with a good upper body control.
- Used for transfers from a sitting or semi-sitting position.
- Used for transfers in mobile and stationary hoists.

### FEATURES

- The sling is used for transferring user from a wheelchair / chair.
- The sling provides good support under the thigh and the back of pelvis to the shoulder blade.
- The sling ensures a good upright posture that facilitates the correct positioning of the user.
- The sling has a large opening around lower body which makes toileting and personal hygiene and changing of clothes easier.
- Split padded leg straps.
- Medium seat opening.



### TECHNICAL SPECIFICATIONS

SWL: 275 kg

Material: Polyester

Sizes: XS-S-M-L-XL-XXL

### ITEM NUMBERS

Low Back sling - polyester

|            | Item number |
|------------|-------------|
| <b>XS</b>  | 009-09633   |
| <b>S</b>   | 009-09634   |
| <b>M</b>   | 009-09635   |
| <b>L</b>   | 009-09636   |
| <b>XL</b>  | 009-09637   |
| <b>XXL</b> | 009-09638   |

## SLING SIZE GUIDE

The sling has a coloured strap on upper back with size and midline mark, which ensures correct placement of the sling.

Colour code is presented on the size chart.

The image below shows a red colour code, which means that the sling is in size S (Small):



|   | XS    | S     | M     | L     | XL    | XXL   |
|---|-------|-------|-------|-------|-------|-------|
| <b>Y</b><br>(hip measurement)                     | 37 cm | 42 cm | 48 cm | 52 cm | 68 cm | 78 cm |
| <b>X</b><br>(measurement from neck to lower back) | 50 cm | 55 cm | 66 cm | 75 cm | 75 cm | 75 cm |



Please consult the technical manual for further specifications. All figures and measurements are for guidance only and may differ in actual situations. WinnCare Nordic reserves the right to changes without notice. Product photos may include accessories. Please consult your dealer for further information.

## FOR MORE INFORMATION VISIT

[winnCare.com](http://winnCare.com)

E-mail: [info@winnCare.dk](mailto:info@winnCare.dk)

Phone: +45 70273720

WIN  CARE