

## **MORFEA**<sup>®</sup>











## **KEY FEATURES**

- ▶ High Resilient foam (HR)
- Density: ≥ 34 Kg/m³
- Hardness: 3 kPa
- Dimensions: 198 x 88 x 14 cm (77.9 x 34.6 x 5.5 inches)
- Mattress weight: 9.5 Kg
- Covers available: Promust PU and Promust PU HD

## KEY BENEFITS

Meet the needs of persons facing a loss of autonomy, discontinuously bedridden and able to mobilize by themselves.

## $({\sf INDICATIONS})$

HELP IN PREVENTING PRESSURE ULCERS FOR PATIENT AT NONE TO LOW RISK OF PRESSURE ULCER DEVELOPMENT.

BEDRIDDEN FEW HOURS A DAY AND ABLE TO MOBILIZING BY HIMSELF PATIENT WEIGHT: 30 TO 120 KG/ 4.7 TO 18.8 ST

- ▶ Single-piece molded High Resilient (HR) polyurethane foam mattress with 4 cut corners
- ▶ Removable and washable integral cover

	DESCRIPTION	ITEM REFERENCES ACCORDING TO PACKAGING			
		CARDBOX	PALLET BOX	ROLLED AND COMPRESSED (RC)	FIRE STANDARDS
	MORFEA® Promust PU cover (not installed)	VMB14/HIP	VMB14/HIP-BOX (10 units)	VMB14RC1-HIP (1 unit) VMB14RC-HIP (5 units)	EN 597-1 and 2 GPEM D1 90 Classe C
	MORFEA® Promust PU HD cover (not installed)	VMB14/HIPH	VMB14/HIPH-BOX (10 units)	_	EN 597-1 and 2 GPEM D1 90 Classe D

2 years warranty







