



An innovative approach

The new range of



visco-elastic curved cushions

A new range of innovative visco-elastic cushions



Help in reducing the risk of pressure ulcers for seated patients

## ELEGANT and DISCRET



**Benefits:** Reduces the risk of pressure ulcers in the ischiatic and coccyx region  
Better adaptation to the users morphology and their material environment

**Recommended usage:** Reduces the risk of pressure ulcers for patient with a “medium to high” risk without stability problems (without pelvic support) or with stability problems (with pelvic support).

**Limitation of use:** Asymetry of the pelvic region (according to medical device).  
Pneumatic individually deformable cells like the Bi-compartment Kineris cushion would be recommended.

**16 available sizes:** 4 for overweight persons and 2 for persons with important ischiatic bony protuberances.


Cushion is more than 6 cm thick under the ischions (ischiatic region)

# A new look



## EFFICIENCY

Comfort  
Abduction of hips  
*Contoured front edge*



Reduces the risk of coccyx pressure ulcers: *sacrum and coccyx cut-out.*

Reduced risks of the patient slipping forward (*tilted sitting position*)

Easier patient transfers:  
*Rounded edges*

Improved prevention of pressure ulcers in ischiatic regions  
*Minimum thickness of the cushions now 6 cm*



## ERGONOMICS

**Easy indentification of the size of the products**

*Size shown on moulded foam*



**Easy cushion and cover orientation**

*Label on the rear of the cover*



**Handle on the base of the cover**

*Less risks of tearing  
More hygienic*

**Secure**

*No slip base on the cover*



# Dimensions

DESIGNATION	REFERENCES	WIDTH CM	DEPTH CM	HEIGHT CM MINI / MAXI*	PATIENT WEIGHT KG (STONE)
<b>ALOVA<sup>®</sup> Contoured Cushion multi-sizes with 2 covers **</b>	VCAG4340	43	40	6 / 9	30 to 90 Kg (4 to 14 st)
	VCAG4540	45	40	6 / 9	30 to 120 Kg (4 to 18 st)
	VCAG3841	38	41	6 / 9	30 to 90 Kg (4 to 14 st)
	VCAG4041	40	41	6 / 9	30 to 90 Kg (4 to 14 st)
	VCAG4341	43	41	6 / 9	30 to 100 Kg (4 to 15 st)
	VCAG4043	40	43	6 / 9	30 to 90 Kg (4 to 14 st)
	VCAG4544	45	44	6 / 9	30 to 120 Kg (4 to 18 st)
	VCAG5044	50	44	6 / 9	30 to 135 Kg (4 to 21 st)
<b>ALOVA<sup>®</sup> Contoured Cushion multi-sizes XL with 2 covers **</b>	VCAG4343XL-E	43	43	9 / 12	30 to 110 (4 to 16 st)
	VCAG5040XL-E	50	40	9 / 12	30 to 135 Kg (4 to 21 st)
<b>ALOVA<sup>®</sup> Contoured Cushion multi-sizes XXL with 2 covers **</b>	VCAG5250	52	50	7 / 10	120 to 160 (18 to 25 st)
	VCAG5450	54	50	7 / 10	120 to 160 (18 to 25 st)
	VCAG5650	56	50	7 / 10	120 to 200 Kg (18 to 31 st)
	VCAG6050	60	50	7 / 10	120 to 200 Kg (18 to 31 st)
<b>ALOVA<sup>®</sup> Contoured with Pelvic support with 2 covers</b>	VCAG3841BP	38	41	6 / 9	30 to 90 Kg (4 to 14 st)
	VCAG4041BP	40	41	6 / 9	30 to 90 Kg (4 to 14 st)
	VCAG4341BP	43	41	6 / 9	30 to 100 Kg (4 to 15 st)
	VCAG4343BP	43	43	6 / 9	30 to 100 Kg (4 to 15 st)
	VCAG4544BP	45	44	6 / 9	30 to 120 Kg (4 to 18 st)
	VCAG5044BP	50	44	6 / 9	30 to 135 Kg (4 to 21 st)
	VCAG4347BP	43	47	6 / 9	30 to 120 Kg (4 to 18 st)

WINNOCARE  
g r o u p