



# TRITON DUAL

Powered, height adjustable, tilt-in-space shower commode chair

## BETTER WORKING ENVIRONMENT AND MORE COMFORT

- Correct working posture eliminates strain and stress on the caregiver's back, shoulders and knees.
- With the powered seat tilt and adjustable backrest, Triton Dual makes it easy for the user to find the most comfortable sitting position.
- The lift height of 40 cm helps both the caregiver and the user.
- The height adjustment allows an optimal positioning of the chair, ensuring a correct working posture for the care giver.
- For the user, the height adjustment means increased independence; it is possible to adjust the chair, so the user can reach the sink and look in the mirror - thereby the user can help take care of his own personal hygiene.

## FEATURES

- Safety and washability (IPX6W rated).
- The digital control box ensures, that if the seat is raised while tilted forwards, the Dual Security function will simultaneously tilt the seat backwards. This way, the user will not have the feeling of falling out of the chair while it is moving.
- Transport function: with one press on the hand control the chair can be put in a position used for safe transport of the user.
- The chair also comes with two memory functions which makes it quick and easy for the carer to adjust the chair in a desired position.

## TECHNICAL SPECIFICATION

Max. user weight: 150 kg  
(330 lbs/23,62 stone)

Distance btw. arm rest: 470 mm

Seat depth: 430 mm

Seat width: 460 mm

Seat height: 535-950 mm

Total width: 570 mm

Length: 880-1050 mm

Total height: 1065-1480 mm (without head rest)

Weight: 41 kg

Tilt angle: +7° til -30° +2°

IP class: IPx6W

Charging time: 6 hours

Battery capacity: 80-100 lifting cycles

**ITEM NUMBER: 044-00050**

## PRODUCT HIGHLIGHTS

- Powered height adjustment (Soft Start/Stop).
- Powered tilt (+7° / -30°) and adjustable backrest for more comfort and user stability.
- IPx6W rated: Can be cleaned in washing machines up to 250 times due to LINAK IPX6 Washable DURA™ (50 washes is standard).
- Developed from the Easy-Clean principle, which means there are no small corners or edges. The chair has closed screws/joinings and fully welded steel, which prevents water from entering the pipes. It is also treated with double coating which makes it rust resistant.
- Dual Security Function (it can not be raised and tilted forward simultaneously)
- Transport function via the hand control for safe transportation of the user.
- Two programmable memory buttons on the hand control for automatic positioning and adjustment.
- Multi adjustable headrest with integrated rod which does not scratch or damage the wall behind the toilet
- Fits over all standard toilets and suitable for use in the shower
- Max user weight 150 kg
- 



Please consult the technical manual for further specifications. All figures and measurements are for guidance only and may differ in actual situations. Winnocare Nordic reserves the right to changes without notice. Product photos may include accessories. Please consult your dealer for further information.

## FOR MORE INFORMATION VISIT

[winnocare.com](http://winnocare.com)

E-mail: [info@winnocare.dk](mailto:info@winnocare.dk)

Telefon: +45 70273720

WIN CARE