

ELEV UP[®]



ELEV UP COMPACT SIT-TO-STAND LIFT

FOR TRANSFERS IN THE SITTING POSITION

- ▶ Ideal for tight spaces, overall width: 62.5 cm.
- ▶ New non-slip foot rest which is very easy to clean.
- ▶ Height of tibia support can be adjusted.
- ▶ Easy to dismantle and transport.
- ▶ Ergonomic gripping handles.




Emergency stop button

SLING FOR THE ELEV UP



Polypropylene sling, washable at 30°

Maintenance:



SIT-TO-STAND LIFT

COMPACT

- ▶ Easy-to-handle and convenient in small spaces.

ERGONOMIC

- ▶ For transfers in the sitting position.

| DESCRIPTION | FOOT LENGTH | OVERALL WIDTH | SPACING | WEIGHT | | MAX PATIENT WEIGHT | CASTORS | STANDARDS | MOTORIZATION |
|---------------------|-------------|---------------|----------------|-----------|----|--------------------|---------|-------------------------------------|-----------------------|
| ELEV UP | 110 cm | 62.5 cm | 50.5 - 78.5 cm | 40 kg | | 150 kg | Ø 75 mm | NF EN 10535 EN 601-1 EN 601-2 | LINAK IP 45 and IP 65 |
| | A | B | C | D | E | | | | |
| POLYPROPYLENE SLING | 90 cm | 35 cm | Max 60 cm | Max 85 kg | 40 | 150 kg | | | |

Remote control, electric and mechanical parts: 5-year warranty
Batteries, slings and weighing unit: 1 year

Photos not contractually binding

WINNOCARE Group - 4, Le Pas du Château - 85670 Saint-Paul-Mont-Penit - France

Tél : +33 (0)2 51 98 55 64 - Fax : +33 (0)2 51 98 59 07 - E-mail : info@medicatlantique.fr - www.winnocare.fr

In the interest of constantly improving its products, the WINNOCARE Group reserves the rights to modify the technical specifications of any product without prior notice

Class I medical device in accordance with European Directive 93/42/EEC

