



WALKING SLING

For users who need extra support and safety during rehabilitation exercises

USAGE

- For users who need support while gait or balance training.
- For training with WinnCare's ErgoTrainer or stationary hoists.
- Strap system on the sling can be applied to most 2-point spreader bars.

FEATURES

- The sling is allowing many customising options.
- It is possible to adjust the spread around the abdomen, the length of the leg straps and shoulder straps.
- The sling ensures safety and support, in case the user lost balance under training.
- Removable padding on the shoulder and leg straps.
- Elastic back and sides, ensuring a good fit for the user.
- Coloured strap on upper back with size and midline mark, which ensure correct placement of the sling.
- Possible to place and buckle it from behind. It is also still possible to buckle it on the front.

TECHNICAL SPECIFICATIONS

SWL: 180 kg

Material: Polyester

Sizes: XXS-XS-S-M-L-XL-XL

ITEM NUMBERS

Walking sling - polyester

	Item number
XXS	009-09572
XS	009-09883
S	009-09884
M	009-09885
L	009-09886
XL	009-09887



SLING SIZE GUIDE

The sling has a coloured strap placed under the label with size indication and midline mark, which ensures correct placement of the sling.

Colour code is presented on the size chart.

The image below shows a yellow colour code, which means that the sling is in size M (Medium):



	XXS	XS	S	M	L	XL
Belly dimension	62 – 75 cm	76 – 84 cm	85 – 92 cm	93 – 105 cm	106 - 114 cm	116 - 130 cm



Please consult the technical manual for further specifications. All figures and measurements are for guidance only and may differ in actual situations. WinnCare Nordic reserves the right to changes without notice. Product photos may include accessories. Please consult your dealer for further information.

FOR MORE INFORMATION VISIT

winnCare.com

E-mail: info@winnCare.dk

Phone: +45 70273720

WIN  CARE